

For Research Use Only

IL-3R alpha/CD123 Protein, Human, Recombinant (His)



Catalog Number: CM19311

产品信息

Catalog Number:
CM19311

CAS号:

分子式:

主要靶点:

分子量:

34.5 kDa (predicted); 50-60 kDa
(reducing condition, due to
glycosylation)

CAS No.:

储存

Shipping with blue ice.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.