

For Research Use Only

PUMA BH3 acetate



www.ptgcn.com

Catalog Number: CM19552

产品信息

Catalog Number:
CM19552

分子量:
3109.35

CAS No.:

CAS号:

分子式:

$C_{130}H_{206}N_{42}O_{45}S$

主要靶点:

BCL

主要通路:

凋亡

靶点活性

p53:26 nM(Kd)

描述

PUMA BH3 acetate is a p53 positive regulator of apoptosis (PUMA) BH3 domain polypeptide that acts as a direct activator of Bak with a Kd of 26 nM.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.