

For Research Use Only

# PerCP Anti-Mouse CD4 (GK1.5)

Catalog Number: CP-65104



## Basic Information

Catalog Number:

CP-65104

Size:

100ug, 200  $\mu$ g/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC039137

GeneID (NCBI):

12504

UNIPROT ID:

P06332

Full Name:

CD4 antigen

Purification Method:

Affinity purification

CloneNo.:

GK1.5

Excitation/Emission maxima  
wavelengths:

482 nm / 677 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Background Information

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.25 ug PerCP Anti-Mouse CD4 (GK1.5) (CP-65104, Clone: GK1.5) or unstained. Cells were then stained with 0.25 ug CoraLite® Plus 647 Anti-Mouse CD3 (17A2) (CL647-65077, Clone: 17A2). Cells were not fixed.