

For Research Use Only

# FITC Plus Anti-Rat CD43 (W3/13) Rabbit Recombinant Antibody



Catalog Number: FITC-65556

## Basic Information

Catalog Number:

FITC-65556

Size:

100ug, 500  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Affinity purification

CloneNo.:

W3/13

Excitation/Emission maxima  
wavelengths:

495 nm / 524 nm

## Applications

Tested Applications:

FC

Species Specificity:

Rat

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> rat splenocytes were surface stained with 0.25 ug FITC Plus Anti-Rat CD43 (W3/13) Rabbit Recombinant Antibody (FITC-65556, Clone: W3/13), or FITC-conjugated Rabbit IgG control Rabbit PolyAb (FITC-30000). Cells were not fixed.