

For Research Use Only

MultiPro™ 5CFLX Anti-Human IL-3RA (Polyclonal)



Catalog Number: G13655-1-5C

Basic Information

Catalog Number:

G13655-1-5C

Size:

500ug/mL

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4452

GenBank Accession Number:

BC035407

GeneID (NCBI):

3563

ENSEMBL Gene ID:

ENSG00000185291

UNIPROT ID:

P26951

Full Name:

MultiPro™ 5CFLX Anti-Human IL-3RA
(Polyclonal)

Conjugate:

5CFLX

Full Oligo Sequence:

CGGAGATGTGTATAAGAGACAGTTCT
CTGGCCAGGCTCCCATATAAGAAA

Barcode Sequence:

TTCTCTGGCCAGGCT

Applications

Tested Applications:

Single Cell (Intra)

Species Specificity:

Human

Background Information

Storage

Storage:

2-8°C

Storage Buffer:

PBS with 1mM EDTA and 0.09% sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



G13655-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G13655-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.