

For Research Use Only

MultiPro™ 5CFLX Anti-Human LEF1 (Polyclonal)



Catalog Number: G14972-1-5C

Basic Information

Catalog Number: G14972-1-5C	GenBank Accession Number: BC050632	Conjugate: 5CFLX
Size: 500ug/mL	GeneID (NCBI): 51176	Full Oligo Sequence: CGGAGATGTGTATAAGAGACAGCATA CAGGCTGACAACCCATATAAGAAA
Source: Rabbit	ENSEMBL Gene ID: ENSG00000138795	Barcode Sequence: CATACAGGCTGACAA
Isotype: IgG	UNIPROT ID: Q9UJU2	
Immunogen Catalog Number: AG6882	Full Name: MultiPro™ 5CFLX Anti-Human LEF1 (Polyclonal)	

Applications

Tested Applications:
Single Cell (Intra)

Species Specificity:
Human

Background Information

Lymphoid enhancer-binding factor 1(LEF1) belongs to a family of regulatory protein share homology with high mobility group protein-1, and it's a nuclear protein expressed in pre-B and T cells. LEF1 has a role in the Wnt signaling pathway and hair cell differentiation and follicle morphogenesis. Together with CTNNB1 and EP300, LEF1 activates transcription of target genes. Isoform 5 transcriptionally activates the fibronectin promoter, binds to and represses transcription from the E-cadherin promoter in a CTNNB1-independent manner, and is involved in reducing cellular aggregation and increasing cell migration of pancreatic cancer cells. Isoform 1 transcriptionally activates MYC and CCND1 expression and enhances proliferation of pancreatic tumor cells. MECs can give rise to seven cell types of the SAE and SMGs following severe airway injury. MECs progressively adopted a basal cell phenotype on the SAE and established lasting progenitors capable of further regeneration following reinjury. MECs activate Wnt-regulated transcription factors (Lef-1/TCF7) following injury and Lef-1 induction in cultured MECs promoted transition to a basal cell phenotype. Surprisingly, dose-dependent MEC conditional activation of Lef-1 in vivo promoted self-limited airway regeneration in the absence of injury. Thus, modulating the Lef-1 transcriptional program in MEC-derived progenitors may have regenerative medicine applications for lung diseases. (<https://doi.org/10.1016/j.stem.2018.03.017>) The phosphorylation may affects LEF1 protein's theoretical molecular weight when tested. 40-70 kD bands have also been reported (PMID: 22261717; 17063141).

Storage

Storage:
2-8°C
Storage Buffer:
PBS with 1mM EDTA and 0.09% sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



G14972-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G14972-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.