

For Research Use Only

MultiPro™ 5CFLX Anti-Human Ezrin (Polyclonal)



Catalog Number: G26056-1-5C

Basic Information

Catalog Number: G26056-1-5C	GenBank Accession Number: BC013903	Conjugate: 5CFLX
Size: 500ug/mL	GeneID (NCBI): 7430	Full Oligo Sequence: CGGAGATGTGTATAAGAGACAGGCCA CCAGTGGCTCCCCATATAAGAAA
Source: Rabbit	ENSEMBL Gene ID: ENSG00000092820	Barcode Sequence: GCCACCAGTGGCTCC
Isotype: IgG	UNIPROT ID: P15311	
Immunogen Catalog Number: AG23530	Full Name: MultiPro™ 5CFLX Anti-Human Ezrin (Polyclonal)	

Applications

Tested Applications:
Single Cell (Intra)

Species Specificity:
Human

Background Information

Storage


Storage:
2-8°C

Storage Buffer:
PBS with 1mM EDTA and 0.09% sodium azide


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



G26056-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G26056-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.