

For Research Use Only

# MultiPro™ 5CFLX Anti-Human CD83 (HB15e)



Catalog Number: G65171-1-5C

## Basic Information

Catalog Number:

G65171-1-5C

Size:

500ug/mL

Source:

Mouse

Isotype:

Mouse IgG1 kappa

GenBank Accession Number:

BC030830

GeneID (NCBI):

9308

ENSEMBL Gene ID:

ENSG00000112149

UNIPROT ID:

Q01151

Full Name:

MultiPro™ 5CFLX Anti-Human CD83  
(HB15e)

CloneNo.:

HB15e

Conjugate:

5CFLX

Full Oligo Sequence:

CGGAGATGTGTATAAGACAGCGGC

CAACTACTGCCCCATATAAGAAA

Barcode Sequence:

CGGCCAACTACTCTG

## Applications

Tested Applications:

Single Cell (Intra)

Species Specificity:

Human

## Background Information

### Storage

Storage:

2-8°C

Storage Buffer:

PBS with 1mM EDTA and 0.09% sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



G65171-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G65171-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.