

For Research Use Only

MultiPro™ 5CFLX Anti-Human TGM2 (2D4C11)



Catalog Number: G68006-1-5C

Basic Information

Catalog Number: G68006-1-5C	GenBank Accession Number: BC003551	CloneNo.: 2D4C11
Size: 500ug/mL	GeneID (NCBI): 7052	Conjugate: 5CFLX
Source: Mouse	ENSEMBL Gene ID: ENSG00000198959	Full Oligo Sequence: CGGAGATGTGTATAAGACAGTATA CTGAAGGATATCCCATATAAGAAA
Isotype: IgG2a	UNIPROT ID: P21980	Barcode Sequence: TATACTGAAGGATAT
Immunogen Catalog Number: AG7462	Full Name: MultiPro™ 5CFLX Anti-Human TGM2 (2D4C11)	

Applications

Tested Applications:
Single Cell (Intra)
Species Specificity:
Human

Background Information

Storage

Storage:
2-8°C
Storage Buffer:
PBS with 1mM EDTA and 0.09% sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



G68006-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G68006-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.