

For Research Use Only

# HRP-conjugated TNFSF18 Polyclonal antibody



Catalog Number: HRP-23899

## Basic Information

Catalog Number:

HRP-23899

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20979

GenBank Accession Number:

BC069319

GeneID (NCBI):

8995

UNIPROT ID:

Q9UNG2

Full Name:

tumor necrosis factor (ligand)  
superfamily, member 18

Calculated MW:

177 aa, 20 kDa

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB

Species Specificity:

human

Positive Controls:

WB : COLO 320 cells, SGC-7901 cells

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with HRP-23899 (TNFSF18 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.