

IHC*easy* AGXT Ready-To-Use IHC Kit

Catalog Number: **KHC0507**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

AGXT (AGT) is a key enzyme in glyoxylate detoxification. Mutations of AGXT lead to primary hyperoxaluria type I, a rare metabolic disease which is characterized by deficient AGXT, excessive hepatic oxalate production that leads to renal failure.


Synonyms

AGT, AGT1, AGXT, AGXT1, PH1, SPAT, SPT, TLH6

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using KHC0507 (AGXT IHC Kit).



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using KHC0507 (AGXT IHC Kit).