

IHC*easy* CLTC Ready-To-Use IHC Kit

Catalog Number: **KHC0687**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Mouse Monoclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

Clathrin is the major protein of the polyhedral coat of coated pits and vesicles which entrap specific macromolecules during receptor-mediated endocytosis. The clathrin molecule has a triskelion shape. Each clathrin triskelion is composed of three identical heavy chains (CLTC) and three light chains of two types, LCA (CLTA) and LCB (CLTB). The heavy chain provides the structural backbone of the clathrin lattice.


Synonyms

CHC, CHC17, Clathrin heavy chain 1, clathrin, heavy chain (Hc), CLH 17, CLH17, CLTC, CLTCL2, Hc, KIAA0034


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using KHC0687 (CLTC IHC Kit).




Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using KHC0687 (CLTC IHC Kit).




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using KHC0687 (CLTC IHC Kit).



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC0687 (CLTC IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using KHC0687 (CLTC IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat colon tissue slide using KHC0687 (CLTC IHC Kit).