

IHC*easy* PSMD11 Ready-To-Use IHC Kit

Catalog Number: **KHC1419**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Rabbit Polyclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

The 26S proteasome is a 2.5-MDa molecular machine that degrades ubiquitinated proteins in eukaryotic cells. It consists of a proteolytic core particle and two 19S regulatory particles (RPs) composed of 6 ATPase (RPT) and 13 non-ATPase (RPN) subunits. PSMD11 gene encodes 19S proteasome subunit RPN6. Increased levels of PSMD11 and a corresponding increased assembly of the 26S/30S proteasome is correlated with high proteasome activity. In vitro ectopic expression of PSMD11 is sufficient to increase proteasome assembly and activity. Latest research has implicated PSMD11 in regulation of human embryonic stem cells function.


Synonyms

p44.5, PSMD11, Rpn6, S9


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC1419 (PSMD11 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using KHC1419 (PSMD11 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human lung tissue slide using KHC1419 (PSMD11 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using KHC1419 (PSMD11 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using KHC1419 (PSMD11 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat lung tissue slide using KHC1419 (PSMD11 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using KHC1419 (PSMD11 IHC Kit).