

IHC*easy* TRIM15 Ready-To-Use IHC Kit

Catalog Number: **KHC1914**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

| Component | Size | Concentration |
|--------------------------|--------------------|---------------|
| Antigen Retrieval Buffer | 100 mL | 50× |
| Washing Buffer | 100 mL ×2 | 20× |
| Blocking Buffer | 5 mL | RTU |
| Primary Antibody | 5 mL | RTU |
| Secondary Antibody | 5 mL | RTU |
| Chromogen Component A | 0.2 mL | RTU |
| Chromogen Component B | 4 mL | RTU |
| Signal Enhancer | 5 mL | RTU |
| Counter Staining Reagent | 5 mL | RTU |
| Mounting Media | 5 mL | RTU |
| Control Slide | 1 slide (Optional) | FFPE |
| Datasheet | 1 Copy | |
| Manual | 1 Copy | |

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

TRIM15 (Tripartite motif-containing protein 15) is also named as RNF93, ZNF178, ZNFB7 and belongs to the TRIM/RBCC family. It acts as an ubiquitin-protein ligase and contributes to the endogenous restriction of retroviruses in cells.


Synonyms

RNF93, tripartite motif containing 15, TRIM, RING finger protein 93, EC:2.3.2.27


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC1914 (TRIM15 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC1914 (TRIM15 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC1914 (TRIM15 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat small intestine tissue slide using KHC1914 (TRIM15 IHC Kit).