

IHC*easy* BHLHE41 Ready-To-Use IHC Kit

Catalog Number: **KHC1926**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

BHLHE41 (basic helix-loop-helix family member e41), also called SHARP1 and DEC2, is a member of the BHLH family proteins that are involved in differentiation, circadian rhythms and the response to hypoxia. In esophageal, osteosarcomas, prostate, breast, and endometrial cancer cells, BHLHE41 functions as an anti-apoptotic factor and promotes cell proliferation and invasiveness. It suppresses cell proliferation and metastasis, and also facilitates pancreatic as well as gastric cancer cell apoptosis. BHLHE41 has an important role in the proliferation of glioblastoma and could serve as a novel candidate for targeted therapy of glioblastoma.


Synonyms

BHLHB3, BHLHE41, DEC2, hDEC2, SHARP 1, SHARP1


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse heart tissue slide using KHC1926 (BHLHE41 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using KHC1926 (BHLHE41 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using KHC1926 (BHLHE41 IHC Kit).