

ADAMDEC1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00007-3

Capture Antibody Information

Catalog Number:
82985-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230367B10
Reactivity:
human
GenBank:
BC074877
Immunogen Catalog Number:
Ag12265

Conjugate:
Unconjugated
Full name:
ADAM-like, decysin 1
Gene ID:
27299

Detection Antibody Information

Catalog Number:
82985-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230367D12
Reactivity:
human
GenBank:
BC074877
Immunogen Catalog Number:
Ag12265

Conjugate:
Unconjugated
Full name:
ADAM-like, decysin 1
Gene ID:
27299

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00007-3 targets ADAMDEC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ADAMDEC1 Recombinant antibody, PBS Only (Capture) 82985-2-PBS (230367B10). 100 µg. Concentration 1 mg/ml.

Detection antibody: ADAMDEC1 Recombinant antibody, PBS Only (Capture/Detector) 82985-5-PBS (230367D12). 100 µg. Concentration 1 mg/ml.

Alternative ADAMDEC1 matched antibody pairs: MP00007-1, MP00007-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00007-3, ADAMDEC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82985-2-PBS. Detection antibody: 82985-5-PBS. Standard: Ag12265. Range: 0.625-20 ng/mL.