

FLT4/VEGFR3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00009-4

Capture Antibody Information

Catalog Number:

85076-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242430C2

Reactivity:

human

GenBank:

NM_182925

Conjugate:

Unconjugated

Full name:

fms-related tyrosine kinase 4

Gene ID:

2324

Detection Antibody Information

Catalog Number:

85076-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242430F3

Reactivity:

human

GenBank:

NM_182925

Conjugate:

Unconjugated

Full name:

fms-related tyrosine kinase 4

Gene ID:

2324

Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00009-4 targets FLT4/VEGFR3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FLT4 Recombinant antibody, PBS Only (Capture) 85076-1-PBS (242430C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: FLT4 Recombinant antibody, PBS Only (Detector) 85076-3-PBS (242430F3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00009-4, FLT4/VEGFR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85076-1-PBS. Detection antibody: 85076-3-PBS. Standard: Eg2494. Range: 0.781-100 ng/mL.