

For Research Use Only

ICAM-1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00028-1

Capture Antibody Information

Catalog Number:
83069-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230342F1

Reactivity:
Human

GenBank:
BC015969

Conjugate:
Unconjugated

Full name:
intercellular adhesion molecule 1

Gene ID:
3383

Detection Antibody Information

Catalog Number:
83069-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230342F9

Reactivity:
Human

GenBank:
BC015969

Conjugate:
Unconjugated

Full name:
intercellular adhesion molecule 1

Gene ID:
3383

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00028-1 targets ICAM-1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ICAM-1 Recombinant antibody, PBS Only (Capture) 83069-1-PBS (230342F1). 100 µg. Concentration 1 mg/mL.

Detection antibody: ICAM-1 Recombinant antibody, PBS Only (Detector) 83069-3-PBS (230342F9). 100 µg. Concentration 1 mg/mL.

Alternative ICAM-1 matched antibody pairs: MP00028-2, MP00028-3, MP00028-4, MP00028-5, MP50085-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00028-1, ICAM-1 Recombinant Matched Antibody Pair - BSA and Azide free, measured by cytometric bead array. Capture antibody: 83069-1-PBS. Detection antibody: 83069-3-PBS. Standard: Eg0570. Range: 0.625-20 ng/mL