

For Research Use Only

# CAMP Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00042-2

## Capture Antibody Information

**Catalog Number:**  
82940-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230157E9  
**Reactivity:**  
Human  
**GenBank:**  
BC055089  
**Immunogen Catalog Number:**  
Ag2622

**Conjugate:**  
Unconjugated  
**Full name:**  
cathelicidin antimicrobial peptide  
**Gene ID:**  
820

## Detection Antibody Information

**Catalog Number:**  
82940-6-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230157H2  
**Reactivity:**  
Human  
**GenBank:**  
BC055089  
**Immunogen Catalog Number:**  
Ag2622

**Conjugate:**  
Unconjugated  
**Full name:**  
cathelicidin antimicrobial peptide  
**Gene ID:**  
820

## Applications

**Tested Applications:**  
Sandwich ELISA

**Range:**  
0.39-25 ng/mL (Sandwich ELISA)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00042-2 targets CAMP in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: CAMP Recombinant antibody, PBS Only (Capture) 82940-4-PBS (230157E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: CAMP Recombinant antibody, PBS Only (Detector) 82940-6-PBS (230157H2). 100 µg. Concentration 1 mg/mL.

Alternative CAMP matched antibody pairs: MP00042-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP00042-2, CAMP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82940-4-PBS. Detection antibody: 82940-6-PBS. Standard: Ag2622. Range: 0.39-25 ng/mL.