

RAB32 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00043-1

Capture Antibody Information

Catalog Number:
82960-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230315F1
Reactivity:
human
GenBank:
BC015061
Immunogen Catalog Number:
Ag1450

Conjugate:
Unconjugated
Full name:
RAB32, member RAS oncogene family
Gene ID:
10981

Detection Antibody Information

Catalog Number:
82960-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230315B11
Reactivity:
human
GenBank:
BC015061
Immunogen Catalog Number:
Ag1450

Conjugate:
Unconjugated
Full name:
RAB32, member RAS oncogene family
Gene ID:
10981

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.


Product Information

MP00043-1 targets RAB32 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: RAB32 Recombinant antibody, PBS Only (Capture/Detector) 82960-4-PBS (230315F1). 100 µg. Concentration 1 mg/ml.
Detection antibody: RAB32 Recombinant antibody, PBS Only (Detector) 82960-2-PBS (230315B11). 100 µg. Concentration 1 mg/ml.
Alternative RAB32 matched antibody pairs: MP00043-2, MP00043-3
Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00043-1, RAB32 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82960-4-PBS. Detection antibody: 82960-2-PBS. Standard: Ag1450. Range: 2.5-20 ng/mL.