

# Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00047-2

## Capture Antibody Information

**Catalog Number:**  
83053-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230292B6  
**Reactivity:**  
Human  
**GenBank:**  
BC127667

**Conjugate:**  
Unconjugated  
**Full name:**  
integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)  
**Gene ID:**  
3690

## Detection Antibody Information

**Catalog Number:**  
83053-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230292B2  
**Reactivity:**  
Human  
**GenBank:**  
BC127667

**Conjugate:**  
Unconjugated  
**Full name:**  
integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)  
**Gene ID:**  
3690

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.312-5 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00047-2 targets Integrin beta-3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Integrin beta 3 Recombinant antibody, PBS Only (Capture) 83053-2-PBS (230292B6). 100 µg. Concentration 1 mg/ml.

Detection antibody: Integrin beta 3 Recombinant antibody, PBS Only (Detector) 83053-1-PBS (230292B2). 100 µg. Concentration 1 mg/ml.

Alternative Integrin beta-3 matched antibody pairs: MP00047-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00047-2, Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83053-2-PBS. Detection antibody: 83053-1-PBS. Standard: SY00181. Range: 0.312-5 ng/mL.