

For Research Use Only

# Mouse PCSK9 Recombinant Matched Antibody Pair, PBS Only



Catalog Number: MP00059-3

## Capture Antibody Information

**Catalog Number:**  
83004-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230536F6  
**Reactivity:**  
Mouse  
**GenBank:**  
NM\_153565.2

**Conjugate:**  
Unconjugated  
**Full name:**  
proprotein convertase  
subtilisin/kexin type 9  
**Gene ID:**  
100102

## Detection Antibody Information

**Catalog Number:**  
83004-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230536C6  
**Reactivity:**  
Mouse  
**GenBank:**  
NM\_153565.2

**Conjugate:**  
Unconjugated  
**Full name:**  
proprotein convertase  
subtilisin/kexin type 9  
**Gene ID:**  
100102

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00059-3 targets PCSK9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse PCSK9 Recombinant antibody, PBS Only (Capture) 83004-3-PBS (230536F6). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: Mouse PCSK9 Recombinant antibody, PBS Only (Capture/Detector) 83004-2-PBS (230536C6). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative PCSK9 matched antibody pairs: MP00059-1, MP00059-2, MP50055-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00059-3, Mouse PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83004-3-PBS. Detection antibody: 83004-2-PBS. Standard: Eg0640. Range: 0.625-20 ng/mL.