

For Research Use Only

IL-1R1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00062-1

Capture Antibody Information

Catalog Number:
82990-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230192E12

Reactivity:
human

GenBank:
BC075063

Immunogen Catalog Number:
Ag25998

Conjugate:
Unconjugated

Full name:
interleukin 1 receptor, type I

Gene ID:
3554

Detection Antibody Information

Catalog Number:
82990-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230192C9

Reactivity:
human

GenBank:
BC075063

Immunogen Catalog Number:
Ag25998

Conjugate:
Unconjugated

Full name:
interleukin 1 receptor, type I

Gene ID:
3554

Applications

Tested Applications:
Sandwich ELISA

Range:
0.188-12 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00062-1 targets IL-1R1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: IL-1R1 Recombinant antibody, PBS Only (Capture) 82990-2-PBS (230192E12). 100 μ g. Concentration 1 mg/ml.

Detection antibody: IL-1R1 Recombinant antibody, PBS Only (Detector) 82990-1-PBS (230192C9). 100 μ g. Concentration 1 mg/ml.

Alternative IL-1R1 matched antibody pairs: MP50175-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00062-1, IL-1R1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82990-2-PBS. Detection antibody: 82990-1-PBS. Standard: Ag25998. Range: 0.188-12 ng/mL