

For Research Use Only

EPOR Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00065-1

Capture Antibody Information

Catalog Number:
83056-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230117G3
Reactivity:
Human
GenBank:
NM_000121.4
Immunogen Catalog Number:
Ag33956

Conjugate:
Unconjugated
Full name:
erythropoietin receptor
Gene ID:
2057

Detection Antibody Information

Catalog Number:
83056-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230117F9
Reactivity:
Human
GenBank:
NM_000121.4
Immunogen Catalog Number:
Ag33956

Conjugate:
Unconjugated
Full name:
erythropoietin receptor
Gene ID:
2057

Applications

Tested Applications:
Sandwich ELISA

Range:
0.625-40 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00065-1 targets EPOR in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: EPOR Recombinant antibody, PBS Only (Capture) 83056-4-PBS (230117G3). 100 µg. Concentration 1 mg/ml.

Detection antibody: EPOR Recombinant antibody, PBS Only (Detector) 83056-3-PBS (230117F9). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00065-1, EPOR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83056-4-PBS. Detection antibody: 83056-3-PBS. Standard: Ag33956. Range: 0.625-40 ng/mL