

PICK1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00068-2

Capture Antibody Information

Catalog Number:
82919-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230185A9
Reactivity:
Human
GenBank:
BC017561
Immunogen Catalog Number:
Ag1443

Conjugate:
Unconjugated
Full name:
protein interacting with PRKCA 1
Gene ID:
9463

Detection Antibody Information

Catalog Number:
82919-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230185B8
Reactivity:
Human, mouse, rat, pig
GenBank:
BC017561
Immunogen Catalog Number:
Ag1443

Conjugate:
Unconjugated
Full name:
protein interacting with PRKCA 1
Gene ID:
9463

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00068-2 targets PICK1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PICK1 Recombinant antibody, PBS Only (Capture) 82919-3-PBS (230185A9). 100 µg. Concentration 1 mg/mL.

Detection antibody: PICK1 Recombinant antibody, PBS Only (Detector) 82919-1-PBS (230185B8). 100 µg. Concentration 1 mg/mL.

Alternative PICK1 matched antibody pairs: MP00068-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Standard curve of MP00068-2, PICK1 Recombinant Matched Antibody Pair - BSA and Azide free, measured by cytometric bead array. Capture antibody: 82919-3-PBS. Detection antibody: 82919-1-PBS. Standard: Ag1443. Range: 0.625-20 ng/mL