

For Research Use Only

ABCD1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00070-2

Capture Antibody Information

Catalog Number:
83038-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230226F6
Reactivity:
Human
GenBank:
BC025358
Immunogen Catalog Number:
Ag0453

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family D (ALD), member 1
Gene ID:
215

Detection Antibody Information

Catalog Number:
83038-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230226G11
Reactivity:
Human
GenBank:
BC025358
Immunogen Catalog Number:
Ag0453

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family D (ALD), member 1
Gene ID:
215

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
0.625-20 ng/mL (Cytometric Bead Array)
78-5000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00070-2 targets ABCD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: ABCD1 Recombinant antibody, PBS Only (Capture) 83038-2-PBS (230226F6). 100 µg. Concentration 1 mg/ml.

Detection antibody: ABCD1 Recombinant antibody, PBS Only (Detector) 83038-3-PBS (230226G11). 100 µg. Concentration 1 mg/ml.

Alternative ABCD1 matched antibody pairs: MP00070-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00070-2, ABCD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83038-2-PBS. Detection antibody: 83038-3-PBS. Standard: Ag0453. Range: 78-5000 pg/mL.



Cytometric bead array standard curve of MP00070-2, ABCD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83038-2-PBS. Detection antibody: 83038-3-PBS. Standard: Ag0453. Range: 0.625-20 ng/mL.