

For Research Use Only

ARHGEF19 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00072-1

Capture Antibody Information

Catalog Number:
82966-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230290B5
Reactivity:
Human
GenBank:
BC040640
Immunogen Catalog Number:
Ag11645

Conjugate:
Unconjugated
Full name:
Rho guanine nucleotide exchange factor (GEF) 19
Gene ID:
128272

Detection Antibody Information

Catalog Number:
82966-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230290D1
Reactivity:
Human
GenBank:
BC040640
Immunogen Catalog Number:
Ag11645

Conjugate:
Unconjugated
Full name:
Rho guanine nucleotide exchange factor (GEF) 19
Gene ID:
128272

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00072-1 targets ARHGEF19 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ARHGEF19 Recombinant antibody, PBS Only (Capture) 82966-3-PBS (230290B5). 100 µg. Concentration 1 mg/mL.

Detection antibody: ARHGEF19 Recombinant antibody, PBS Only (Detector) 82966-4-PBS (230290D1). 100 µg. Concentration 1 mg/mL.

Alternative ARHGEF19 matched antibody pairs: MP00072-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00072-1, ARHGEF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82966-3-PBS. Detection antibody: 82966-4-PBS. Standard: Ag11645. Range: 0.625-20 ng/mL.