

EFNA2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00081-2

Capture Antibody Information

Catalog Number:
82951-5-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230251H4

Reactivity:
Human

GenBank:
NM_001405

Conjugate:
Unconjugated

Full name:
ephrin-A2

Gene ID:
1943

Detection Antibody Information

Catalog Number:
82951-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230251D2

Reactivity:
Human

GenBank:
NM_001405

Conjugate:
Unconjugated

Full name:
ephrin-A2

Gene ID:
1943

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00081-2 targets EFNA2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EFNA2 Recombinant antibody, PBS Only (Capture) 82951-5-PBS (230251H4). 100 µg. Concentration 1 mg/mL.

Detection antibody: EFNA2 Recombinant antibody, PBS Only (Detector) 82951-3-PBS (230251D2). 100 µg. Concentration 1 mg/mL.

Alternative EFNA2 matched antibody pairs: MP00081-1, MP00081-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00081-2, EFNA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82951-5-PBS. Detection antibody: 82951-3-PBS. Standard: SY00429. Range: 0.313-10 ng/mL