

ABCA13 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00084-1

Capture Antibody Information

Catalog Number:
83060-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230288F5
Reactivity:
Human
GenBank:
NM_152701

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family A (ABC1), member 13
Gene ID:
154664

Detection Antibody Information

Catalog Number:
83060-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230288G4
Reactivity:
Human
GenBank:
NM_152701

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family A (ABC1), member 13
Gene ID:
154664

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00084-1 targets ABCA13 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ABCA13 Recombinant antibody, PBS Only (Capture) 83060-1-PBS (230288F5). 100 µg. Concentration 1 mg/mL.

Detection antibody: ABCA13 Recombinant antibody, PBS Only (Detector) 83060-2-PBS (230288G4). 100 µg. Concentration 1 mg/mL.

Alternative ABCA13 matched antibody pairs: MP00084-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00084-1, ABCA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83060-1-PBS. Detection antibody: 83060-2-PBS. Standard: SY00509. Range: 0.312-20 ng/mL.