

PAK7 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00086-3

Capture Antibody Information

Catalog Number:
82935-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230201A6
Reactivity:
Human
GenBank:
BC024179
Immunogen Catalog Number:
Ag29190

Conjugate:
Unconjugated
Full name:
p21 protein (Cdc42/Rac)-activated kinase 7
Gene ID:
57144

Detection Antibody Information

Catalog Number:
82935-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230201A2
Reactivity:
Human
GenBank:
BC024179
Immunogen Catalog Number:
Ag29190

Conjugate:
Unconjugated
Full name:
p21 protein (Cdc42/Rac)-activated kinase 7
Gene ID:
57144

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00086-3 targets PAK7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PAK7 Recombinant antibody, PBS Only (Capture) 82935-3-PBS (230201A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: PAK7 Recombinant antibody, PBS Only (Capture/Detector) 82935-2-PBS (230201A2). 100 µg. Concentration 1 mg/mL.

Alternative PAK7 matched antibody pairs: MP00086-1, MP00086-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00086-3, PAK7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82935-3-PBS. Detection antibody: 82935-2-PBS. Standard: Ag29190. Range: 0.313-20 ng/mL