

For Research Use Only

HNRNPC Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00088-1

Capture Antibody Information

Catalog Number:
82798-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230373F6
Reactivity:
Human
GenBank:
BC003394
Immunogen Catalog Number:
Ag2356

Conjugate:
Unconjugated
Full name:
heterogeneous nuclear ribonucleoprotein C (C1/C2)
Gene ID:
3183

Detection Antibody Information

Catalog Number:
82798-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230373D7
Reactivity:
Human
GenBank:
BC003394
Immunogen Catalog Number:
Ag2356

Conjugate:
Unconjugated
Full name:
heterogeneous nuclear ribonucleoprotein C (C1/C2)
Gene ID:
3183

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00088-1 targets HNRNPC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HNRNPC Recombinant antibody, PBS Only (Capture) 82798-5-PBS (230373F6). 100 µg. Concentration 1 mg/ml.

Detection antibody: HNRNPC Recombinant antibody, PBS Only (Detector) 82798-4-PBS (230373D7). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00088-1, HNRNPC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82798-5-PBS. Detection antibody: 82798-4-PBS. Standard: Ag2356. Range: 1.25-20 ng/mL.