

For Research Use Only

ENG Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00089-1

Capture Antibody Information

Catalog Number:
82993-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230047C1

Reactivity:
Human

GenBank:
BC014271

Conjugate:
Unconjugated

Full name:
endoglin

Gene ID:
2022

Detection Antibody Information

Catalog Number:
82993-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230047C5

Reactivity:
Human

GenBank:
BC014271

Conjugate:
Unconjugated

Full name:
endoglin

Gene ID:
2022

Applications

Tested Applications:
Sandwich ELISA

Range:
0.156-10 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00089-1 targets ENG in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Endoglin/CD105 Recombinant antibody, PBS Only (Capture) 82993-1-PBS (230047C1). 100 µg. Concentration 1 mg/ml.

Detection antibody: Endoglin/CD105 Recombinant antibody, PBS Only (Detector) 82993-2-PBS (230047C5). 100 µg. Concentration 1 mg/ml.

Alternative ENG matched antibody pairs: MP00089-2, MP50103-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00089-1, ENG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82993-1-PBS. Detection antibody: 82993-2-PBS. Standard: Eg0587. Range: 0.156-10 ng/mL