

DDX39 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00111-2

Capture Antibody Information

Catalog Number:
83083-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230375B12
Reactivity:
human
GenBank:
BC032128
Immunogen Catalog Number:
Ag2311

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-Asp) box polypeptide 39
Gene ID:
10212

Detection Antibody Information

Catalog Number:
83083-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230375E8
Reactivity:
human
GenBank:
BC032128
Immunogen Catalog Number:
Ag2311

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-Asp) box polypeptide 39
Gene ID:
10212

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00111-2 targets DDX39 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDX39 Recombinant antibody, PBS Only (Capture) 83083-1-PBS (230375B12). 100 µg. Concentration 1 mg/ml.

Detection antibody: DDX39 Recombinant antibody, PBS Only (Detector) 83083-2-PBS (230375E8). 100 µg. Concentration 1 mg/ml.

Alternative DDX39 matched antibody pairs: MP00111-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00111-2, DDX39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83083-1-PBS. Detection antibody: 83083-2-PBS. Standard: Ag2311. Range: 1.25-20 ng/mL