

ABCA6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00118-1

Capture Antibody Information

Catalog Number:
83084-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230223A11
Reactivity:
Human
GenBank:
NM_080284
Immunogen Catalog Number:
Ag29524

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family A (ABC1), member 6
Gene ID:
23460

Detection Antibody Information

Catalog Number:
83084-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230223E4
Reactivity:
Human
GenBank:
NM_080284
Immunogen Catalog Number:
Ag29524

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family A (ABC1), member 6
Gene ID:
23460

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00118-1 targets ABCA6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ABCA6 Recombinant antibody, PBS Only (Capture) 83084-1-PBS (230223A11). 100 µg. Concentration 1 mg/ml.

Detection antibody: ABCA6 Recombinant antibody, PBS Only (Detector) 83084-2-PBS (230223E4). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00118-1, ABCA6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83084-1-PBS. Detection antibody: 83084-2-PBS. Standard: Ag29524. Range: 2.5-20 ng/mL.