

For Research Use Only

TBCD Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00119-1

Capture Antibody Information

Catalog Number: 82911-5-PBS	Clone ID: 230209H2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: tubulin folding cofactor D
Isotype: IgG	GenBank: BC003094	Gene ID: 6904
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6310	

Detection Antibody Information

Catalog Number: 82911-4-PBS	Clone ID: 230209G2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: tubulin folding cofactor D
Isotype: IgG	GenBank: BC003094	Gene ID: 6904
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6310	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	---------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP00119-1 targets TBCD in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TBCD Recombinant antibody, PBS Only (Capture) 82911-5-PBS (230209H2). 100 µg. Concentration 1 mg/ml.

Detection antibody: TBCD Recombinant antibody, PBS Only (Detector) 82911-4-PBS (230209G2). 100 µg. Concentration 1 mg/ml.

Alternative TBCD matched antibody pairs: MP00119-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00119-1, TBCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82911-5-PBS. Detection antibody: 82911-4-PBS. Standard: Ag6310. Range: 0.625-20 ng/mL