

For Research Use Only

TBCD Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00119-2

Capture Antibody Information

Catalog Number:
82911-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230209A2

Reactivity:
Human

GenBank:
BC003094

Immunogen Catalog Number:
Ag6310

Conjugate:
Unconjugated

Full name:
tubulin folding cofactor D

Gene ID:
6904

Detection Antibody Information

Catalog Number:
82911-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230209G2

Reactivity:
Human

GenBank:
BC003094

Immunogen Catalog Number:
Ag6310

Conjugate:
Unconjugated

Full name:
tubulin folding cofactor D

Gene ID:
6904

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00119-2 targets TBCD in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TBCD Recombinant antibody, PBS Only (Capture) 82911-3-PBS (230209A2). 100 µg. Concentration 1 mg/mL.

Detection antibody: TBCD Recombinant antibody, PBS Only (Detector) 82911-4-PBS (230209G2). 100 µg. Concentration 1 mg/mL.

Alternative TBCD matched antibody pairs: MP00119-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00119-2, TBCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82911-3-PBS. Detection antibody: 82911-4-PBS. Standard: Ag6310. Range: 0.625-20 ng/mL