

For Research Use Only

NFX1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00127-1

Capture Antibody Information

Catalog Number:
82924-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230211A10
Reactivity:
Human
GenBank:
BC012151
Immunogen Catalog Number:
Ag9448

Conjugate:
Unconjugated
Full name:
nuclear transcription factor, X-box binding 1
Gene ID:
4799

Detection Antibody Information

Catalog Number:
82924-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230211D1
Reactivity:
Human
GenBank:
BC012151
Immunogen Catalog Number:
Ag9448

Conjugate:
Unconjugated
Full name:
nuclear transcription factor, X-box binding 1
Gene ID:
4799

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00127-1 targets NFX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFX1 Recombinant antibody, PBS Only (Capture) 82924-3-PBS (230211A10). 100 µg. Concentration 1 mg/mL.

Detection antibody: NFX1 Recombinant antibody, PBS Only (Detector) 82924-4-PBS (230211D1). 100 µg. Concentration 1 mg/mL.

Alternative NFX1 matched antibody pairs: MP00127-2, MP00127-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00127-1, NFX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82924-3-PBS. Detection antibody: 82924-4-PBS. Standard: Ag9448. Range: 0.625-20 ng/mL