

SCN3B Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00132-2

Capture Antibody Information

Catalog Number:
82959-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230316A4
Reactivity:
human
GenBank:
BC126265
Immunogen Catalog Number:
Ag23201

Conjugate:
Unconjugated
Full name:
sodium channel, voltage-gated, type III, beta
Gene ID:
55800

Detection Antibody Information

Catalog Number:
82959-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230316H5
Reactivity:
human
GenBank:
BC126265
Immunogen Catalog Number:
Ag23201

Conjugate:
Unconjugated
Full name:
sodium channel, voltage-gated, type III, beta
Gene ID:
55800

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00132-2 targets SCN3B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SCN3B Recombinant antibody, PBS Only (Capture) 82959-2-PBS (230316A4). 100 µg. Concentration 1 mg/ml.

Detection antibody: SCN3B Recombinant antibody, PBS Only (Detector) 82959-5-PBS (230316H5). 100 µg. Concentration 1 mg/ml.

Alternative SCN3B matched antibody pairs: MP00132-1, MP00132-3, MP00132-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00132-2, SCN3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82959-2-PBS. Detection antibody: 82959-5-PBS. Standard: Ag23201. Range: 0.625-20 ng/mL