

INTS5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00141-1

Capture Antibody Information

Catalog Number:
83154-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230462D7
Reactivity:
Human
GenBank:
BC060841
Immunogen Catalog Number:
Ag5129

Conjugate:
Unconjugated
Full name:
integrator complex subunit 5
Gene ID:
80789

Detection Antibody Information

Catalog Number:
83154-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230462F2
Reactivity:
Human
GenBank:
BC060841
Immunogen Catalog Number:
Ag5129

Conjugate:
Unconjugated
Full name:
integrator complex subunit 5
Gene ID:
80789

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00141-1 targets INTS5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: INTS5 Recombinant antibody, PBS Only (Capture) 83154-4-PBS (230462D7). 100 µg. Concentration 1 mg/ml.

Detection antibody: INTS5 Recombinant antibody, PBS Only (Capture/Detector) 83154-3-PBS (230462F2). 100 µg. Concentration 1 mg/ml.

Alternative INTS5 matched antibody pairs: MP00141-2, MP00141-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00141-1, INTS5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83154-4-PBS. Detection antibody: 83154-3-PBS. Standard: Ag5129. Range: 2.5-20 ng/mL.