

For Research Use Only

PGF Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00146-1

Capture Antibody Information

Catalog Number:

83170-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230532H3

Reactivity:

Human

GenBank:

BC007789

Immunogen Catalog Number:

Ag0969

Conjugate:

Unconjugated

Full name:

placental growth factor

Gene ID:

5228

Detection Antibody Information

Catalog Number:

83170-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230532C12

Reactivity:

Human

GenBank:

BC007789

Immunogen Catalog Number:

Ag0969

Conjugate:

Unconjugated

Full name:

placental growth factor

Gene ID:

5228

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00146-1 targets PGF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PGF Recombinant antibody, PBS Only (Capture) 83170-4-PBS (230532H3). 100 μ g. Concentration 1 mg/mL.

Detection antibody: PGF Recombinant antibody, PBS Only (Capture/Detector) 83170-2-PBS (230532C12). 100 μ g. Concentration 1 mg/mL.

Alternative PGF matched antibody pairs: MP00146-2, MP00146-3, MP50063-1, MP50063-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00146-1, PGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83170-4-PBS. Detection antibody: 83170-2-PBS. Standard: Ag0969. Range: 0.625-20 ng/mL