

nectin 4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00147-1

Capture Antibody Information

Catalog Number:

83223-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240009E4

Reactivity:

Human

GenBank:

BC010423

Conjugate:

Unconjugated

Full name:

poliovirus receptor-related 4

Gene ID:

81607

Affinity: $K_D = 8.49 \times 10^{-10}$ $K_{Off} = 2.20 \times 10^{-4}$ $K_{On} = 2.59 \times 10^5$

Detection Antibody Information

Catalog Number:

83223-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240009F2

Reactivity:

Human

GenBank:

BC010423

Conjugate:

Unconjugated

Full name:

poliovirus receptor-related 4

Gene ID:

81607

Affinity: $K_D = 2.28 \times 10^{-10}$ $K_{Off} = 2.98 \times 10^{-5}$ $K_{On} = 2.62 \times 10^5$

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00147-1 targets nectin 4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: nectin 4 Recombinant antibody, PBS Only (Capture) 83223-1-PBS (240009E4). 100 µg. Concentration 1 mg/ml.

Detection antibody: nectin 4 Recombinant antibody, PBS Only (Detector) 83223-4-PBS (240009F2). 100 µg. Concentration 1 mg/ml.

Alternative nectin 4 matched antibody pairs: MP00147-2, MP00147-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.comW: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00147-1, nectin 4 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83223-1-PBS. Detection antibody: 83223-4-PBS. Standard: Eg0415. Range: 0.625-10 ng/mL.