

For Research Use Only

KLF9 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00148-3

Capture Antibody Information

Catalog Number: 83183-2-PBS	Clone ID: 230436G4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: Kruppel-like factor 9
Isotype: IgG	GenBank: BC074880	Gene ID: 687
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34411	

Detection Antibody Information

Catalog Number: 83183-1-PBS	Clone ID: 230436B1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: Kruppel-like factor 9
Isotype: IgG	GenBank: BC074880	Gene ID: 687
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34411	

Applications

Tested Applications: Cytometric bead array	Range: 1.25-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	--------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP00148-3 targets KLF9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KLF9 Recombinant antibody, PBS Only (Capture) 83183-2-PBS (230436G4). 100 µg. Concentration 1 mg/ml.

Detection antibody: KLF9 Recombinant antibody, PBS Only (Detector) 83183-1-PBS (230436B1). 100 µg. Concentration 1 mg/ml.

Alternative KLF9 matched antibody pairs: MP00148-1, MP00148-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00148-3, KLF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83183-2-PBS. Detection antibody: 83183-1-PBS. Standard: Ag34411. Range: 1.25-20 ng/mL