

For Research Use Only

HOXB6 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00150-2

Capture Antibody Information

Catalog Number:
83155-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230433C11
Reactivity:
Human
GenBank:
BC014651
Immunogen Catalog Number:
Ag14189

Conjugate:
Unconjugated
Full name:
homeobox B6
Gene ID:
3216

Detection Antibody Information

Catalog Number:
83155-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230433G1
Reactivity:
Human
GenBank:
BC014651
Immunogen Catalog Number:
Ag14189

Conjugate:
Unconjugated
Full name:
homeobox B6
Gene ID:
3216

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00150-2 targets HOXB6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HOXB6 Recombinant antibody, PBS Only (Capture) 83155-1-PBS (230433C11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: HOXB6 Recombinant antibody, PBS Only (Detector) 83155-3-PBS (230433G1). 100 μ g. Concentration 1 mg/mL.

Alternative HOXB6 matched antibody pairs: MP00150-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00150-2, HOXB6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83155-1-PBS. Detection antibody: 83155-3-PBS. Standard: Ag14189. Range: 2.5-20 ng/mL.