

For Research Use Only

CBX5 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00181-2

Capture Antibody Information

Catalog Number:
83258-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240145E5
Reactivity:
Human
GenBank:
BC006821
Immunogen Catalog Number:
Ag34810

Conjugate:
Unconjugated
Full name:
chromobox homolog 5 (HP1 alpha homolog, Drosophila)
Gene ID:
23468

Detection Antibody Information

Catalog Number:
83258-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240145D4
Reactivity:
Human
GenBank:
BC006821
Immunogen Catalog Number:
Ag34810

Conjugate:
Unconjugated
Full name:
chromobox homolog 5 (HP1 alpha homolog, Drosophila)
Gene ID:
23468

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00181-2 targets CBX5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CBX5 Recombinant antibody, PBS Only (Capture) 83258-3-PBS (240145E5). 100 µg. Concentration 1 mg/ml.

Detection antibody: CBX5 Recombinant antibody, PBS Only (Detector) 83258-1-PBS (240145D4). 100 µg. Concentration 1 mg/ml.

Alternative CBX5 matched antibody pairs: MP00181-1, MP00181-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00181-2, CBX5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83258-3-PBS. Detection antibody: 83258-1-PBS. Standard: Ag34810. Range: 3.125-100 ng/mL.