

For Research Use Only

# Histone H2B Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00183-1

## Capture Antibody Information

Catalog Number:  
83171-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
230539F7

Reactivity:  
Human

GenBank:  
BC005827

Conjugate:  
Unconjugated

Full name:  
histone cluster 2, H2be

Gene ID:  
8349

## Detection Antibody Information

Catalog Number:  
83171-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
230539H2

Reactivity:  
Human

GenBank:  
BC005827

Conjugate:  
Unconjugated

Full name:  
histone cluster 2, H2be

Gene ID:  
8349

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.625-5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00183-1 targets Histone H2B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Histone H2B Recombinant antibody, PBS Only (Capture) 83171-2-PBS (230539F7). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: Histone H2B Recombinant antibody, PBS Only (Detector) 83171-3-PBS (230539H2). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative Histone H2B matched antibody pairs: MP00183-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00183-1, Histone H2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83171-2-PBS. Detection antibody: 83171-3-PBS. Standard: SY02145Ac. Range: 0.625-5 ng/mL.