

For Research Use Only

SULT1B1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00184-2

Capture Antibody Information

Catalog Number:
83228-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240119F6
Reactivity:
Human
GenBank:
BC010895
Immunogen Catalog Number:
Ag9009

Conjugate:
Unconjugated
Full name:
sulfotransferase family, cytosolic, 1B, member 1
Gene ID:
27284

Detection Antibody Information

Catalog Number:
83228-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240119F12
Reactivity:
Human
GenBank:
BC010895
Immunogen Catalog Number:
Ag9009

Conjugate:
Unconjugated
Full name:
sulfotransferase family, cytosolic, 1B, member 1
Gene ID:
27284

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00184-2 targets SULT1B1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: SULT1B1 Recombinant antibody, PBS Only (Capture/Detector) 83228-2-PBS (240119F6). 100 μ g. Concentration 1 mg/ml.
Detection antibody: SULT1B1 Recombinant antibody, PBS Only (Detector) 83228-3-PBS (240119F12). 100 μ g. Concentration 1 mg/ml.
Alternative SULT1B1 matched antibody pairs: MP00184-1
Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00184-2, SULT1B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83228-2-PBS. Detection antibody: 83228-3-PBS. Standard: Ag9009. Range: 2.5-80 ng/mL.