

For Research Use Only

# SP8 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00190-2

## Capture Antibody Information

<b>Catalog Number:</b> 83230-4-PBS	<b>Clone ID:</b> 240011B12	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> Sp8 transcription factor
<b>Isotype:</b> IgG	<b>GenBank:</b> BC038669	<b>Gene ID:</b> 221833
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag22934	

## Detection Antibody Information

<b>Catalog Number:</b> 83230-1-PBS	<b>Clone ID:</b> 240011A1	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> Sp8 transcription factor
<b>Isotype:</b> IgG	<b>GenBank:</b> BC038669	<b>Gene ID:</b> 221833
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag22934	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 2.5-20 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP00190-2 targets SP8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SP8 Recombinant antibody, PBS Only (Capture) 83230-4-PBS (240011B12). 100 µg. Concentration 1 mg/ml.

Detection antibody: SP8 Recombinant antibody, PBS Only (Capture/Detector) 83230-1-PBS (240011A1). 100 µg. Concentration 1 mg/ml.

Alternative SP8 matched antibody pairs: MP00190-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00190-2, SP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83230-4-PBS. Detection antibody: 83230-1-PBS. Standard: Ag22934. Range: 2.5-20 ng/mL.