

For Research Use Only

PDGFB Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00198-2

Capture Antibody Information

Catalog Number:
83214-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230486G11
Reactivity:
Human
GenBank:
BC029822
Immunogen Catalog Number:
Ag25986

Conjugate:
Unconjugated
Full name:
platelet-derived growth factor beta polypeptide (simian sarcoma viral (v-sis) oncogene homolog)
Gene ID:
5155
Affinity:
 $K_D = 1.16 \times 10^{-9}$
 $K_{Off} = 1.97 \times 10^{-4}$
 $K_{On} = 1.70 \times 10^5$

Detection Antibody Information

Catalog Number:
83214-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230486C9
Reactivity:
Human
GenBank:
BC029822
Immunogen Catalog Number:
Ag25986

Conjugate:
Unconjugated
Full name:
platelet-derived growth factor beta polypeptide (simian sarcoma viral (v-sis) oncogene homolog)
Gene ID:
5155

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00198-2 targets PDGFB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDGFB Recombinant antibody, PBS Only (Capture) 83214-4-PBS (230486G11). 100 µg. Concentration 1 mg/ml.

Detection antibody: PDGFB Recombinant antibody, PBS Only (Detector) 83214-3-PBS (230486C9). 100 µg. Concentration 1 mg/ml.

Alternative PDGFB matched antibody pairs: MP00198-1, MP00198-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00198-2, PDGFB Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83214-4-PBS. Detection antibody: 83214-3-PBS. Standard: Ag25986. Range: 0.625-20 ng/mL