

For Research Use Only

# TYK2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00203-1

## Capture Antibody Information

**Catalog Number:**  
83161-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230176F8  
**Reactivity:**  
Human  
**GenBank:**  
BC014243  
**Immunogen Catalog Number:**  
Ag28457

**Conjugate:**  
Unconjugated  
**Full name:**  
tyrosine kinase 2  
**Gene ID:**  
7297

## Detection Antibody Information

**Catalog Number:**  
83161-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230176C9  
**Reactivity:**  
Human  
**GenBank:**  
BC014243  
**Immunogen Catalog Number:**  
Ag28457

**Conjugate:**  
Unconjugated  
**Full name:**  
tyrosine kinase 2  
**Gene ID:**  
7297

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00203-1 targets TYK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TYK2 Recombinant antibody, PBS Only (Capture/Detector) 83161-3-PBS (230176F8). 100 µg. Concentration 1 mg/ml.

Detection antibody: TYK2 Recombinant antibody, PBS Only (Detector) 83161-2-PBS (230176C9). 100 µg. Concentration 1 mg/ml.

Alternative TYK2 matched antibody pairs: MP00203-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00203-1, TYK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83161-3-PBS. Detection antibody: 83161-2-PBS. Standard: Ag28457. Range: 0.625-20 ng/mL