

# NUDCD2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00205-1

## Capture Antibody Information

<b>Catalog Number:</b> 83217-3-PBS	<b>Clone ID:</b> 230513C5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> NudC domain containing 2
<b>Isotype:</b> IgG	<b>GenBank:</b> BC017934	<b>Gene ID:</b> 134492
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag15579	

## Detection Antibody Information

<b>Catalog Number:</b> 83217-1-PBS	<b>Clone ID:</b> 230513A7	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> NudC domain containing 2
<b>Isotype:</b> IgG	<b>GenBank:</b> BC017934	<b>Gene ID:</b> 134492
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag15579	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 1.25-40 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP00205-1 targets NUDCD2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NUDCD2 Recombinant antibody, PBS Only (Capture) 83217-3-PBS (230513C5). 100 µg. Concentration 1 mg/ml.

Detection antibody: NUDCD2 Recombinant antibody, PBS Only (Detector) 83217-1-PBS (230513A7). 100 µg. Concentration 1 mg/ml.

Alternative NUDCD2 matched antibody pairs: MP00205-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00205-1, NUDCD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83217-3-PBS. Detection antibody: 83217-1-PBS. Standard: Ag15579. Range: 1.25-40 ng/mL.